

# GRID CABINET

*Designed for use in Radiographic Tables and Wall stands*

The QUEST Grid Cabinet (Q-1717GC) is designed to fit into most radiographic tables and wallstands that currently house a standard style grid cabinet.

This grid cabinet accepts a variety of X-Ray grids and most ion chambers.



QUEST Grid Cabinet (Q-1717GC)

The QUEST Cassette Tray (Q-1717CT) is a great complement for this Grid Cabinet - Sold Separately.



*Shown with QUEST Cassette Tray (Q-1717CT) & CR Cassette Sold Separately*

## KEY FEATURES

CAT# Q-1717GC

### GRID CABINET

- Tray locks provide for safe placement and removal of detector
- Accepts X-Ray grids (of 17 1/4" x 18 7/8")
- Accepts ion chambers (of 18" x 18 1/2")
- For Radiographic Tables and Wall Stands

*\*Specifications are subject to change*



**Provides Accurate Positioning For  
Portrait & Landscape Alignment**

[www.questxray.com](http://www.questxray.com)

