


**“DESIGNED TO FIT” YOUR “EXISTING X-RAY TABLE”  
TO PROVIDE SAFER TABLE EXAMINATION CONDITIONS**



Contact QUEST for compatibility to existing X-Ray tables, as some x-ray table designs may require modifications for attachment.

NOTE: QUEST'S "Table Safety Rails" (Q-TSR) are designed to provide additional safety for a patient when on the x-ray table. It is the responsibility of the operator to always assure that all proper safety measures are implemented for each patient.

## OPTIONAL: SAFETY RAILS WALL BRACKET

### SAFETY RAILS WALL BRACKET (Q-TSR-W)

- A wall mounted bracket designed for holding a set of two table safety rails when not in use.

**OPTIMAL FOR ALL "OEM's"  
PROVIDING EXTRA VALUE TO CUSTOMERS**

## KEY FEATURES

CAT# Q-TSR

### TABLE SAFETY RAILS

- Safety Rail: Set of 2 detachable protector devices.
- Height above table = 4.5" (Overall 6" H x 22" L)
- Designed for providing additional patient safety.
- Mounts to front and back of the X-Ray table, Using the X-Ray tables existing support rails.
- Fast and flexible positioning along the entire length of the X-Ray table rail.
- Attaches by sliding onto the existing table side-channel opening and then securing with a simple hand lock.
- Easy to attach and remove, as needed.
- Durable, Lightweight and easy to clean.

*\*Specifications are subject to change*

[www.questxray.com](http://www.questxray.com)



**Provides Safer Table  
Examination Conditions**